


ESys

Complete MRI Compatible In-Room Entertainment System



ESys

Complete MRI Compatible In-Room Entertainment System provides high definition In-magnet room audio and visual presentation for the enjoyment of the patient during an MR exam. May be positioned at either end of the bore. Patient views screen using coil mirrors. Image may be inverted as needed. System complete with stereo headphones and conversation microphone.



 **Invivo**

- High Definition In-Magnet Room Display
- 2-Way Audio, Visual Presentation System
- Plays Video, Audio, MP3's, Other Medias
- MR Safe, RF Shielded, 30" LCD Display
- Compatible at 1.5T and 3.0T
- Easy Roll, Adjustable, Counter Balanced Stand

©2008 Invivo Corporation All Rights Reserved
580207 Rev 01

12501 Research Parkway
Orlando, FL 32826 USA
tel 407-275-3220
tel 800-331-3220
fax 407-249-2022

info@invivocorp.com
www.invivocorp.com